



Blood Cell Storage, Inc.
 4203 23rd Ave W.
 Seattle WA, 98199 USA

Certificate of Compliance

Date of Manufacture: 25 February 2022
 Product Name: BCSI pH SAFE Reader (pH1000™)
 Part Number: 400001
 Lot Number: 2202049†
 BCSI pH1000 Serial Number: BCS0809P0112

BCSI pH1000 QC Standard Serial Number: B396

Serial Number	Fr	Fp	Rp
B396	1966	3194	1.578



B3961966319415780254

NOTE: These criteria are specific for the pH1000 paired with the QC Standard as stated above. The Fr and Fp values must be within 5% of the listed value. The Rp value must be within 1% of the listed value.

USER: _____ EXPIRY†: _____
 Date QC Standard implemented (barcode scanned) One (1) year from date of barcode scan

USER: Record the expiry date on the SQC label.

BCSI certifies that the above referenced product:

- Was tested by BCSI whose Quality Management System is approved by an accredited registering body to the appropriate ISO 13485:2016 Medical Devices Quality Systems Standard and ISO 9001:2015 Quality Systems Standard.
- Was tested in accordance with applicable BCSI Standard Procedures.
- Was tested and complies with visual and performance specifications as specified in the product specification.



Steve Geelhood
 Business Director

25 Feb 22
 Date

400001-07