



Blood Cell Storage, Inc.
 4203 23rd Ave W.
 Seattle WA, 98199 USA

Certificate of Compliance

Date of Manufacture: 19 May 2022

Product Name: BCSI pH SAFE Reader (pH1000™)

Part Number: 400001

Lot Number: 2205022†

BCSI pH1000 Serial Number: BCS0804P0090

BCSI pH1000 QC Standard Serial Number: B440

Serial Number	Fr	Fp	Rp
B440	2215	2539	1.213



B4402215253912130241

NOTE: These criteria are specific for the pH1000 paired with the QC Standard as stated above. The Fr and Fp values must be within 5% of the listed value. The Rp value must be within 1% of the listed value.

USER: _____ EXPIRY†: _____
 Date QC Standard implemented (barcode scanned) One (1) year from date of barcode scan

USER: **Record the expiry date on the SQC label.**

BCSI certifies that the above referenced product:

- Was tested by BCSI whose Quality Management System is approved by an accredited registering body to the appropriate ISO 13485:2016 Medical Devices Quality Systems Standard and ISO 9001:2015 Quality Systems Standard.
- Was tested in accordance with applicable BCSI Standard Procedures.
- Was tested and complies with visual and performance specifications as specified in the product specification.



Steve Geelhood
 Business Director

19 May 2022
 Date

400001-07